

Serum free Cell Cryopreservation Solution 5% DMSO

Product Description:

This newly developed cryopreservation medium is a chemically defined cell preservation solution. All raw materials are **USP-grade** (compliant with US Pharmacopeia standards). Notably, it contains only 5% DMSO.

Application:

This product is specialized for diverse cell lines, with validated performance in primary cells including (but not limited to):

- 293T, A549, K562, U937
- Jurkat, NK,DC,B Lymphocytes
- Peripheral blood mononuclear cells (PBMCs)

Features:

Ready-to-use | Stable storage | Easy operation

Chemically defined | Animal/Human origin-free | Protein-free

No programmed cooling | Simplified freezing protocol

High cell recovery >90% | Ensures experimental reproducibility

Long-term cryostability | Compatible with -80°C & LN₂ storage



Product information:

F106-20: 20ml F106-100: 100ml

Shelf Life: 3 years

Storage Conditions:

2-8°C, light protection recommended.

Instruction for use:

- 1. Collect exponentially growing cells (determined by growth curve analysis)
- 2. Wash the cells once with 1xPBS or the corresponding cell wash solution, and remove the supernatant by centrifugation.
- 3. Add an appropriate amount of cell cryopreservation solution according to the number of cells, and slowly, gently resuspend the cells, and transfer to a cell cryopreservation tube.

Cell Type	Number of Cells/ml
Cell line	1-5*10 ⁶
PBMC	1-5*10 ⁷
Mesenchymal stem cells	1-5*106

- 4. Place the cell cryovial in an appropriate freezing container (eg, a cell cryopreservation box containing isopropanol is better for cryopreservation), and place it at 4 degrees or room temperature for 10-15 minutes (this step is for Make sure that the DMSO has enough time to enter the cell).
- 5. Store in the -80 degree refrigerator.
- Reminder: If the cells need to be frozen in a liquid nitrogen tank, the cells should be stored in a -80 degree refrigerator for at least 24 hours.